



**ORTHO
SOLUTIONS**
IMPLANTATE TECHNOLOGIE



VET - Line





**ORTHO
SOLUTIONS**
IMPLANTATE TECHNOLOGIE

Ortho Solutions GmbH - Eltstraße 2 - D-78573 Wurmlingen
Phone +49(0)7461 - 96357 – 0 - Fax +49(0)7461 - 96357 - 99
www.orthosolutions.de - info@orthosolutions.de



Inhalt / Content

1,5MM MINI-FRAGMENT PLATTE	8
1,5MM MINI-PLATE, STRAIGHT	8
1,5MM T-PLATTE	8
1,5MM T-PLATE	8
1,5 MM H-PLATTE	8
1,5 MM H-PLATE	8
2,0 MM MINI DCP-PLATTE, GERADE	9
2,0 MM MINI DCP PLATE, STRAIGHT	9
2,0 MM KOMPRESSIONS-PLATTE, GERADE	9
2,0 MM COMPRESSION PLATE, STRAIGHT	9
2,0 MM L-PLATTE, 90° LINKS	9
2,0 MM L-PLATE, 90° LEFT	9
2,0 MM L-PLATTE, 90° RECHTS	10
2,0 MM L-PLATE, 90° RIGHT	10
2,0 MM T-PLATTE	10
2,0 MM T-PLATE.....	10
2,0 MM L-PLATTE, SCHRÄG LINKS	10
2,0 MM L-PLATE OBLIQUE LEFT	10
2,0 MM L-PLATTE, SCHRÄG RECHTS	10
2,0 MM L-PLATE OBLIQUE RIGHT.....	10
2,0 MM T-PLATTE	11
2,0 MM T-PLATE.....	11
2,0 MM H-PLATTE	11
2,0 MM H-PLATE	11
2,0 MM MINI DCP-PLATTE	11
2,0 MM MINI DCP PLATE.....	11
2,0 MM MINI DCP-PLATTE	12
2,0 MM MINI DCP PLATE.....	12
2,0 MM ADAPTIONS-PLATTE	12
2,0 MM ADAPTION-PLATE	12
2,0 MM LC/DCP MINI PLATTE	12
2,4 MM DCP MINI PLATTE	13
2,4 MM DCP MINI PLATE.....	13
2,7 MM DCP KOMPRESSIONS-PLATTE	14
2,7 MM SMALL FRAGMENT DCP PLATE.....	14



2,0/2,7 MM ZUSCHNEIDEPLATTE	14
2,0/2,7 MM CUTTING PLATE	14
2,7 MM L-PLATTE, 90° LINKS	15
2,7 MM L-PLATE, 90° LEFT	15
2,7 MM L-PLATTE, 90° RECHTS	15
2,7 MM L-PLATE, 90° RIGHT	15
2,7 MM T-PLATTE	15
2,7 MM T-PLATE.....	15
2,7 MM L-PLATTE, SCHRÄG LINKS	16
2,7 MM L-PLATE OBLIQUE LEFT	16
2,7 MM L-PLATTE, SCHRÄG RECHTS	16
2,7 MM L-PLATE OBLIQUE RIGHT.....	16
2,7 MM MULTIFRAGMENT-PLATTE	16
2,7 MM MULTI FRAGMENT-PLATE	16
2,7 MM LC/DCP KOMPRESSIONS-PLATTE	17
2,7 MM LC/DCP COMPRESSION PLATE.....	17
3,5 MM DCP KOMPRESSIONS-PLATTE	18
3,5 MM DCP COMPRESSION PLATE	18
3,5 MM LC/DCP KOMPRESSIONS-PLATTE	18
3,5 MM LC/DCP COMPRESSION PLATE.....	18
3,5 MM DRITTELROHR PLATTE	19
3,5 MM 1/3 TUBULAR PLATE	19
3,5 MM DRITTELROHR PLATTE MIT KRAGEN	19
3,5 MM 1/3 TUBULAR PLATE WITH COLLAR.....	19
3,5 MM VERSTÄRKTE KLEINFRAGMENT PLATTE	20
3,5 MM STRONG SMALL FRAGMENT PLATE.....	20
3,5 MM VERSTÄRKTE LC-KLEINFRAGMENT PLATTE	20
3,5 MM STRONG LC SMALL FRAGMENT PLATE	20
4,5 MM DCP KOMPRESSIONS-PLATTE, SCHMAL	21
4,5 MM DCP COMPRESSION PLATE SMALL.....	21
4,5 MM LC/DCP KOMPRESSIONS-PLATTE, SCHMAL	22
4,5 MM LC/DCP COMPRESSION PLATE SMALL	22
4,5 MM LC/DCP KOMPRESSIONS-PLATTE, BREIT	22
4,5 MM LC/DCP COMPRESSION PLATE BROAD.....	22
4,5 MM; DCP KOMPRESSIONS-PLATTE, BREIT	23
4,5 MM; DCP COMPRESSION PLATE BROAD.....	23
4,5 MM GROÙE T-PLATTE	23



4,5 MM T-PLATE LARGE	23
2,7 MM REKONSTRUKTIONS PLATTE, GEBOGEN.....	24
2,7 MM RECONSTRUCTIVE PLATE, CURVED.....	24
2,7 MM REKONSTRUKTIONS PLATTE, GERADE	24
2,7 MM RECONSTRUCTIVE PLATE, STRAIGHT.....	24
3,5 MM REKONSTRUKTIONS PLATTE, GERADE	25
3,5 MM RECONSTRUCTIVE PLATE, STRAIGHT.....	25
4,5 MM REKONSTRUKTIONS PLATTE, GERADE	25
4,5 M RECONSTRUCTIVE PLATE, STRAIGHT	25
2,0 MM VETERINÄR PLATTE LINKS.....	26
2,0 MM VETERINARY PLATE LEFT.....	26
2,0 MM VETERINÄR PLATTE RECHTS.....	26
2,0 MM VETERINARY PLATE RIGHT.....	26
1,5 / 2,0 MM VETERINÄR T-PLATTE	26
1,5 / 2,0 MM VETERINARY T-PLATE	26
2,0 MM ACETABULUM PLATTE, GEBOGEN	26
2,0 MM ACETABULUM PLATE, CURVED.....	26
2,0 MM MINI DCP VERLÄNGERUNGSPLATTE	27
2,0 MM MINI DCP LENGTHENING PLATE.....	27
2,4/2,7 MM TPLO-PLATTE.....	27
2,4/2,7 MM TPLO PLATE	27
2,0 MM MINI DCP VERLÄNGERUNGSPLATTE	27
2,0 MM MINI DCP LENGHTENING PLATE.....	27
2,7 MM VETERINÄR PLATTE	28
2,7 MM VETERINÄR PLATE	28
2,0 / 2,7 MM VETERINÄR T-PLATTEN.....	28
2,0 / 2,7 MM VETERINÄR T-PLATES	28
2,7 / 3,5 MM VETERINÄR T-PLATTEN.....	28
2,7 / 3,5 MM VETERINÄR T-PLATES	28
2,7 MM ACETABULUM PLATTE, GEBOGEN	29
2,7 MM ACETABULUM PLATE, CURVED.....	29
3,5 MM MPL-VETERINÄR PLATTE RECHTS	29
3,5 MM MPL-VETERINARY PLATE RIGHT	29
3,5 MM DCP KLEINFRAGMENTPLATTE.....	29
3,5 MM DCP SMALL FRAGMENT PLATE.....	29
3,5 MM VET -ABSTÜTZPLATTE PLATTE.....	30
3,5 MM BUTRESS PLATE	30



3,5 MM VET -ABSTÜTZPLATTE PLATTE.....	30
3,5 MM BUTRESS PLATE	30
3,5 MM TPLO-PLATTE	30
3,5 MM TPLO PLATE.....	30
2,7 MM TPLO-PLATTE	30
2,7 MM TPLO PLATE.....	30
3,5 MM TPLO-PLATTE	31
3,5 MM TPLO PLATE.....	31
3,5 MM TPLO-PLATTE	31
3,5 MM TPLO PLATE.....	31
4,5 MM TPLO-PLATTE	31
4,5 MM TPLO PLATE.....	31
1,5 / 2,0 MM ARTHRODESE PLATTE	31
1,5 / 2,0 MM ARTHRODESIS PLATE	31
2,0 / 2,0 MM ARTHRODESE PLATTE	32
2,0 / 2,0 MM ARTHRODESIS PLATE	32
2,0 / 2,7 MM ARTHRODESE PLATTE	32
2,0 / 2,7 MM ARTHRODESIS PLATE	32
2,0 / 2,7 MM ARTHRODESE PLATTE	32
2,0 / 2,7 MM ARTHRODESIS PLATE	32
2,7 / 3,5 MM ARTHRODESE PLATTE	32
2,7 / 3,5 MM ARTHRODESIS PLATE	32
2,7 / 3,5 MM ARTHRODESE PLATTE	32
2,7 / 3,5 MM ARTHRODESIS PLATE	32
2,7 / 3,5 MM; ARTHRODESE PLATTE.....	33
2,7 / 3,5 MM; ARTHRODESIS PLATE	33
2,7 / 3,5 MM; ARTHRODESE PLATTE.....	33
2,7 / 3,5 MM; ARTHRODESIS PLATE	33
3,5 / 3,5 MM; ARTHRODESE PLATTE.....	33
3,5 / 3,5 MM; ARTHRODESIS PLATE	33
3,5 / 3,5 MM ARTHRODESE PLATTE	33
3,5 / 3,5 MM ARTHRODESIS PLATE	33
1,5 MM KORTIKALIS SCHRAUBE, INNENSECHSKANT.....	34
1,5 MM CORTICAL SCREW, HEXAGONAL RECESS	34
1,5 MM KORTIKALIS SCHRAUBE, KREUZSCHLITZ.....	34
1,5 MM CORTICAL SCREW, CRUCIFORM RECESS.....	34
2,0 MM KORTIKALIS SCHRAUBE, INNENSECHSKANT.....	35



2,0 MM CORTICAL SCREW, HEXAGONAL RECESS	35
2,0 MM KORTIKALIS SCHRAUBE, KREUZSCHLITZ.....	36
2,0 MM CORTICAL SCREW, CRUCIFORM RECESS	36
2,4 MM KORTIKALIS SCHRAUBE, INNENSECHSKANT.....	37
2,4 MM CORTICAL SCREW, HEXAGONAL RECESS	37
2,7 MM KORTIKALIS SCHRAUBE, INNENSECHSKANT.....	38
2,7 MM CORTICAL SCREW, HEXAGONAL RECESS	38
3,5 MM KORTIKALIS SCHRAUBE, INNENSECHSKANT.....	39
3,5 MM CORTICAL SCREW, HEXAGONAL RECESS	39
4,5 MM KORTIKALIS SCHRAUBE, INNENSECHSKANT.....	40
4,5 MM CORTICAL SCREW, HEXAGONAL RECESS	40
2,7 MM VERRIEGELUNGSSCHRAUBE, TORX.....	42
2,7 MM LOCKING SCREW, STARDRIVE	42
2,7 MM VERRIEGELUNGSSCHRAUBE, TORX, POLYAXIAL.....	43
2,7 MM LOCKING SCREW, STARDRIVE, POLYAXIAL	43
3,5 MM VERRIEGELUNGSSCHRAUBE, TORX.....	44
3,5 MM LOCKING SCREW, STARDRIVE	44
3,5 MM VERRIEGELUNGSSCHRAUBE, TORX, POLYAXIAL.....	45
3,5 MM LOCKING SCREW, STARDRIVE, POLYAXIAL	45
3,5 MM SPONGIOSA SCHRAUBEN, INNENSECHSKANT.....	46
3,5 MM CANCELLOUS SCREWS, HEXAGONAL RECESS	46
4,0 MM SPONGIOSA SCHRAUBEN, INNENSECHSKANT.....	47
4,0 MM CANCELLOUS SCREWS, HEXAGONAL RECESS	47



Mini Fragment Platten / mini fragment plates

1,5mm Mini-Fragment Platte

1,5mm Mini-Plate, straight

Material (material): ISO 5832-1
Plattendicke (thickness): 1,0 mm
Plattenbreite (width): 3,8 mm
Lochabstand (hole spacing): 5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6010-0003	15	3
6010-0004	20	4
6010-0005	25	5
6010-0006	30	6
6010-0020	100	20



1,5mm T-Platte

1,5mm T-Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 1,0 mm
Plattenbreite (width): 3,8 mm
Lochabstand (hole spacing): 5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher (Schaft) holes (shaft)	Löcher (Kopf) holes (head)
6010-0339	50	9	3
6010-0349	50	9	4



1,5 mm H-Platte

1,5 mm H-Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 1,0 mm
Plattenbreite (width): 8,2 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6010-0804	9	4





Mini Fragment Platten / mini fragment plates

2,0 mm Mini DCP-Platte, gerade

2,0 mm Mini DCP Plate, straight

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0003	17	3
6011-0004	23	4
6011-0005	29	5
6011-0006	35	6
6011-0007	41	7
6011-0008	47	8
6011-0009	53	9
6011-0020	119	20
6011-0048	251	50



2,0 mm Kompressions-Platte, gerade

2,0 mm Compression Plate, straight

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0022	100	20



2,0 mm L-Platte, 90° links

2,0 mm L-Plate, 90° left

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0102	18	2



Mini Fragment Platten / mini fragment plates

2,0 mm L-Platte, 90° rechts

2,0 mm L-Plate, 90° right

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0202	18	2



2,0 mm T-Platte

2,0 mm T-Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0302	18	2



2,0 mm L-Platte, schräg links

2,0 mm L-Plate oblique left

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0402	20	2



2,0 mm L-Platte, schräg rechts

2,0 mm L-Plate oblique right

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0502	20	2





Mini Fragment Platten / mini fragment plates

2,0 mm T-Platte

2,0 mm T-Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	5,0 mm
Lochabstand (hole spacing):	5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher (Schaft) holes (shaft)	Löcher (Kopf) holes (head)
6011-0338	45	8	3
6011-0339	50	9	3
6011-0349	50	9	4



2,0 mm H-Platte

2,0 mm H-Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	11,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-0804	9	4



2,0 mm Mini DCP-Platte

2,0 mm Mini DCP Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,0 mm
Plattenbreite (width):	5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-1004	22	4
6011-1005	27	5
6011-1006	32	6
6011-1007	37	7
6011-1008	42	8



Mini Fragment Platten / mini fragment plates

2,0 mm Mini DCP-Platte

2,0 mm Mini DCP Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 1,5 mm
 Plattenbreite (width): 5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-1014	22	4
6011-1015	27	5
6011-1016	32	6
6011-1017	37	7
6011-1018	42	8
6011-1019	47	9
6011-1020	52	10
6011-1021	57	11
6011-1022	62	12
6011-1030	102	30



2,0 mm Adaptions-Platte

2,0 mm Adaption-Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 1,0 mm
 Plattenbreite (width): 5,2 mm
 Lochabstand (hole spacing): 5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-1420	100	20



2,0 mm LC/DCP Mini Platte

2,0 mm LC/DCP Mini Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 1,0 mm
 Plattenbreite (width): 4,0 mm
 Lochabstand (hole spacing): 6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-2004	23	4
6011-2005	29	5
6011-2006	35	6
6011-2007	41	7
6011-2008	47	8
6011-2009	53	9
6011-2010	59	10
6011-2020	119	20





2,4 mm DCP Mini Platte

2,4 mm DCP Mini Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 2,3 mm
Plattenbreite (width): 6,5 mm
Lochabstand (hole spacing): 7,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6011-6004	31	4
6011-6005	38	5
6011-6006	45	6
6011-6007	52	7
6011-6008	59	8
6011-6009	66	9
6011-6010	73	10
6011-6011	80	11
6011-6012	87	12
6011-6014	101	14





2,7 mm DCP Kompressions-Platte

2,7 mm Small Fragment DCP Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	
2-6 Loch (hole)	2,0 mm
7-12-Loch (hole)	2,5 mm
Plattenbreite (width):	8,0 mm
Lochabstand (hole spacing):	8,0 mm
Lochabstand Mitte (hole spacing middle):	11,0 mm
bei (for) 6012-0002-L	16,8 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0003	28	3
6012-0004	36	4
6012-0005	44	5
6012-0006	52	6
6012-0007	60	7
6012-0008	68	8
6012-0009	76	9
6012-0010	84	10
6012-0011	92	11
6012-0012	100	12
6012-0013	108	13
6012-0014	116	14



2,0/2,7 mm Zuschneideplatte

2,0/2,7 mm Cutting Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,5 mm
Plattenbreite (width):	7,3 mm
Lochabstand (hole spacing):	6,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0050	300	50





Klein Fragment Platten / small fragment plates

2,7 mm L-Platte, 90° links

2,7 mm L-Plate, 90° left

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	7,0 mm
Lochabstand (hole spacing):	8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0103	32	3



2,7 mm L-Platte, 90° rechts

2,7 mm L-Plate, 90° right

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	7,0 mm
Lochabstand (hole spacing):	8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0203	32	3



2,7 mm T-Platte

2,7 mm T-Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	7,0 mm
Lochabstand (hole spacing):	8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0303	32	3
6012-0304	40	4
6012-0305	48	5
6012-0306	56	6
6012-0307	64	7
6012-0308	72	8



Klein Fragment Platten / small fragment plates

2,7 mm L-Platte, schräg links

2,7 mm L-Plate oblique left

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	7,0 mm
Lochabstand (hole spacing):	8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0403	34,5	3



2,7 mm L-Platte, schräg rechts

2,7 mm L-Plate oblique right

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	7,0 mm
Lochabstand (hole spacing):	8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0503	34,5	3



2,7 mm Multifragment-Platte

2,7 mm Multi Fragment-Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 mm
Plattenbreite (width):	14,5 mm
Lochabstand (hole spacing):	8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-0907	44,5	7





2,7 mm LC/DCP Kompressions-Platte

2,7 mm LC/DCP Compression Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	
2-6 Loch (hole)	2,0 mm
7-12-Loch (hole)	2,5 mm
Plattenbreite (width):	8,0 mm
Lochabstand (hole spacing):	8,0 mm
Lochabstand Mitte (hole spacing middle):	11,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6012-2002	20	2
6012-2003	28	3
6012-2004	36	4
6012-2005	44	5
6012-2006	52	6
6012-2007	60	7
6012-2008	68	8
6012-2009	76	9
6012-2010	84	10
6012-2011	92	11
6012-2012	100	12
6012-2014	116	14



Klein Fragment Platten / small fragment plates

3,5 mm DCP Kompressions-Platte

3,5 mm DCP Compression Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,4 mm
Plattenbreite (width):	10,0 mm
Lochabstand (hole spacing):	12,0 mm
Lochabstand Mitte (hole spacing middle):	16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6020-0002	26	2
6020-0003	38	3
6020-0004	50	4
6020-0005	62	5
6020-0006	74	6
6020-0007	86	7
6020-0008	98	8
6020-0009	110	9
6020-0010	122	10
6020-0011	134	11
6020-0012	146	12
6020-0013	158	13
6020-0014	170	14
6020-0016	194	16
6020-0018	218	18
6020-0020	242	20
6020-0022	266	22



3,5 mm LC/DCP Kompressions-Platte

3,5 mm LC/DCP Compression Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	4,5 mm
Plattenbreite (width):	11,5 mm
Lochabstand (hole spacing):	13,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6020-2004	51	4
6020-2005	64	5
6020-2006	77	6
6020-2007	90	7
6020-2008	103	8
6020-2009	116	9
6020-2010	129	10
6020-2011	142	11
6020-2012	155	12





Klein Fragment Platten / small fragment plates

3,5 mm Drittelrohr Platte

3,5 mm 1/3 Tubular Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,1 mm
Plattenbreite (width):	9,0 mm
Lochabstand (hole spacing):	12,0 mm
Lochabstand Mitte (hole spacing middle):	16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6022-0002	25	2
6022-0003	37	3
6022-0004	49	4
6022-0005	61	5
6022-0006	73	6
6022-0007	85	7
6022-0008	97	8
6022-0009	109	9
6022-0010	121	10
6022-0011	133	11
6022-0012	145	12
6022-0013	157	13
6022-0014	169	14
6022-0015	181	15
6022-0016	193	16
6022-0018	217	18



3,5 mm Drittelrohr Platte mit Kragen

3,5 mm 1/3 Tubular Plate with collar

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,1 mm
Plattenbreite (width):	9,0 mm
Lochabstand (hole spacing):	12,0 mm
Lochabstand Mitte (hole spacing middle):	16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6022-1002	26	2
6022-1003	38	3
6022-1004	50	4
6022-1005	62	5
6022-1006	74	6
6022-1007	86	7
6022-1008	98	8
6022-1009	110	9
6022-1010	122	10
6022-1011	134	11
6022-1012	146	12



Klein Fragment Platten / small fragment plates

3,5 mm verstärkte Kleinfragment Platte

3,5 mm strong small Fragment Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,6 mm
Plattenbreite (width):	12,0 mm
Lochabstand (hole spacing):	12,0 mm
Lochabstand Mitte (hole spacing middle):	16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6024-0005	62	5
6024-0006	74	6
6024-0007	86	7
6024-0008	98	8
6024-0009	110	9
6024-0010	122	10
6024-0011	134	11
6024-0012	146	12
6024-0013	158	13
6024-0014	170	14
6024-0016	194	15
6024-0018	218	18
6024-0020	242	20
6024-0022	266	22



3,5 mm verstärkte LC-Kleinfragment Platte

3,5 mm strong LC small Fragment Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,6 mm
Plattenbreite (width):	12,0 mm
Lochabstand (hole spacing):	12,0 mm
Lochabstand Mitte (hole spacing middle):	13,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6024-2007	86	7
6024-2008	98	8
6024-2009	110	9
6024-2010	122	10
6024-2011	134	11
6024-2012	146	12
6024-2014	170	14
6024-2016	194	16



4,5 mm DCP Kompressions-Platte, schmal

4,5 mm DCP Compression Plate small

Material (material):	ISO 5832-1
Plattendicke (thickness):	4,0 mm
Plattenbreite (width):	12,0 mm
Lochabstand (hole spacing):	16,0 mm
Lochabstand Mitte (hole spacing middle):	25,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6030-0002	39	2
6030-0003	55	3
6030-0004	71	4
6030-0005	87	5
6030-0006	103	6
6030-0007	119	7
6030-0008	135	8
6030-0009	151	9
6030-0010	167	10
6030-0011	183	11
6030-0012	199	12
6030-0013	215	13
6030-0014	231	14
6030-0015	247	15
6030-0016	263	16
6030-0018	295	18
6030-0020	327	20
6030-0022	359	22
6030-0024	391	24



Groß Fragment Platten / large fragment plates

4,5 mm LC/DCP Kompressions-Platte, schmal

4,5 mm LC/DCP Compression Plate small

Material (material):	ISO 5832-1
Plattendicke (thickness):	4,2 mm
Plattenbreite (width):	13,5 mm
Lochabstand (hole spacing):	18,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6030-2002	34	2
6030-2003	52	3
6030-2004	70	4
6030-2005	88	5
6030-2006	106	6
6030-2007	124	7
6030-2008	142	8
6030-2009	160	9
6030-2010	178	10
6030-2011	196	11
6030-2012	214	12
6030-2013	232	13
6030-2014	250	14
6030-2015	268	15
6030-2016	286	16



4,5 mm LC/DCP Kompressions-Platte, breit

4,5 mm LC/DCP Compression Plate broad

Material (material):	ISO 5832-1
Plattendicke (thickness):	5,2 mm
Plattenbreite (width):	17,5 mm
Lochabstand (hole spacing):	18,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6031-2006	106	6
6031-2007	124	7
6031-2008	142	8
6031-2009	160	9
6031-2010	178	10
6031-2011	196	11
6031-2012	214	12
6031-2013	232	13
6031-2014	250	14
6031-2015	268	15
6031-2016	286	16
6031-2017	304	17
6031-2018	322	18





Groß Fragment Platten / large fragment plates

4,5 mm; DCP Kompressions-Platte, breit

4,5 mm; DCP Compression Plate broad

Material (material):	ISO 5832-1
Plattendicke (thickness):	5,0 mm
Plattenbreite (width):	16,2 mm
Lochabstand (hole spacing):	16,0 mm
Lochabstand Mitte (hole spacing middle):	25,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6031-0002	39	2
6031-0003	55	3
6031-0004	71	4
6031-0005	87	5
6031-0006	103	6
6031-0007	119	7
6031-0008	135	8
6031-0009	151	9
6031-0010	167	10
6031-0011	183	11
6031-0012	199	12
6031-0013	215	13
6031-0014	231	14
6031-0015	247	15
6031-0016	263	16
6031-0017	279	17
6031-0018	295	18
6031-0020	327	20
6031-0022	359	22
6031-0024	391	24
6031-0026	423	26



4,5 mm große T-Platte

4,5 mm T-Plate large

Material (material):	ISO 5832-1
Plattendicke (thickness):	2,0 mm
Plattenbreite (width):	17,0 mm
Lochabstand (hole spacing):	33,0 mm
Lochabstand Mitte (hole spacing middle):	16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6034-0303	68	3
6034-0304	84	4
6034-0305	100	5
6034-0306	116	6
6034-0307	134	7
6034-0308	148	8
6034-0309	164	9
6034-0310	180	10
6034-0311	196	11
6034-0312	212	12





2,7 mm Rekonstruktions Platte, gebogen

2,7 mm Reconstructive Plate, curved

Material (material): ISO 5832-1
Plattendicke (thickness): 3,1 mm
Plattenbreite (width): 8,0 mm
Lochabstand (hole spacing): 8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6050-0106	70	5
6050-0108	94	6
6050-0110	118	8
6050-0112	142	10
6050-0114	166	12
6050-0116	190	14
6050-0118	214	16



2,7 mm Rekonstruktions Platte, gerade

2,7 mm Reconstructive Plate, straight

Material (material): ISO 5832-1
Plattendicke (thickness): 3,1 mm
Plattenbreite (width): 8,0 mm
Lochabstand (hole spacing): 8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6051-0005	40	5
6051-0006	48	6
6051-0008	64	8
6051-0010	80	10
6051-0012	96	12
6051-0014	112	14
6051-0016	128	16
6051-0018	144	18
6051-0020	160	20
6051-0022	176	22
6051-0024	192	24



3,5 mm Rekonstruktions Platte, gerade

3,5 mm Reconstructive Plate, straight

Material (material): ISO 5832-1
 Plattendicke (thickness): 2,8 mm
 Plattenbreite (width): 10,0 mm
 Lochabstand (hole spacing): 12,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6050-0003	34	3
6050-0004	46	4
6050-0005	58	5
6050-0006	70	6
6050-0007	82	7
6050-0008	94	8
6050-0009	106	9
6050-0010	118	10
6050-0011	130	11
6050-0012	142	12
6050-0013	154	13
6050-0014	166	14
6050-0015	178	15
6050-0016	190	16
6050-0017	202	17
6050-0018	214	18
6050-0020	238	20
6050-0022	262	22



4,5 mm Rekonstruktions Platte, gerade

4,5 m Reconstructive Plate, straight

Material (material): ISO 5832-1
 Plattendicke (thickness): 2,9 mm
 Plattenbreite (width): 12,0 mm
 Lochabstand (hole spacing): 16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
6052-0006	93	6
6052-0008	125	8
6052-0010	157	10
6052-0012	189	12
6052-0014	221	14
6052-0016	253	16





2,0 mm Veterinär Platte links

2,0 mm Veterinary Plate left

Material (material): ISO 5832-1
Plattendicke (thickness): 1,5 mm
Plattenbreite (width): 5,0 mm
Kopfbreite (head width): 10,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8011-0103	28	3



2,0 mm Veterinär Platte rechts

2,0 mm Veterinary Plate right

Material (material): ISO 5832-1
Plattendicke (thickness): 1,5 mm
Plattenbreite (width): 5,0 mm
Kopfbreite (head width): 10,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8011-0203	28	3



1,5 / 2,0 mm Veterinär T-Platte

1,5 / 2,0 mm Veterinary T-Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 1,5 mm
Plattenbreite (width): 5,0 mm
Kopfbreite (head width): 11,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8011-0302	20	2
8011-0304	32	4
8011-0305	38	5
8011-0306	44	6
8011-0307	50	7
8011-0308	56	8



2,0 mm Acetabulum Platte, gebogen

2,0 mm Acetabulum Plate, curved

Material (material): ISO 5832-1
Plattendicke (thickness): 1,5 mm
Plattenbreite (width): 5,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8011-0606	32	6





2,0 mm Mini DCP Verlängerungsplatte

2,0 mm Mini DCP Lengthening Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 1,5 mm
Plattenbreite (width): 5,3 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8011-5004	42	4
8011-5006	52	6
8011-5008	62	8
8011-5010	72	10
8011-5012	82	12
8011-5014	92	14
8011-5016	102	16



2,4/2,7 mm TPLO-Platte

2,4/2,7 mm TPLO Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 2,0 mm
Plattenbreite (width): 7,2 mm
Kopfbreite (head width): 13,5 mm

Art.Nr. Item no.	Länge mm length mm	Lochabstand mm hole spacing mm	Löcher Holes
8011-6306	34	52	3



2,0 mm Mini DCP Verlängerungsplatte

2,0 mm Mini DCP Lengthening Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 1,5 mm
Plattenbreite (width): 5,3 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8011-5137	92	32
8011-5147	102	42
8011-5157	112	52



Veterinär Platten / veterinary plates

2,7 mm Veterinär Platte

2,7 mm Veterinär Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	2,5 mm
Plattenbreite (width):	8,0 mm
Kopfbreite (head width):	16,4 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes	Ausrichtung direction
8012-0103	41	3	rechts /right
8012-0203	41	3	links / left



2,0 / 2,7 mm Veterinär T-Platten

2,0 / 2,7 mm Veterinär T-Plates

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,5 mm
Plattenbreite (width):	8,0 mm
Kopfbreite (head width):	12,4 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8012-0304	36	4
8012-0305	44	5
8012-0306	52	6



2,7 / 3,5 mm Veterinär T-Platten

2,7 / 3,5 mm Veterinär T-Plates

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,0 mm
Plattenbreite (width):	9,4 mm
Kopfbreite (head width):	18,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8012-0505	65	5
8012-0506	75	6
8012-0507	85	7
8012-0509	105	9





2,7 mm Acetabulum Platte, gebogen

2,7 mm Acetabulum Plate, curved

Material (material): ISO 5832-1
Plattendicke (thickness): 2,0 mm
Plattenbreite (width): 8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8012-0604	39	4
8012-0605	47	5
8012-0606	53	6



3,5 mm MPL-Veterinär Platte rechts

3,5 mm MPL-Veterinary Plate right

Material (material): ISO 5832-1
Plattendicke (thickness): 4,3 mm
Plattenbreite (width): 10,0 mm
Kopfbreite (head width): 17,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes	Ausrichtung direction
8020-0103	64	3	rechts / right
8020-0203	64	3	links /left



3,5 mm DCP Kleinfragmentplatte

3,5 mm DCP small fragment plate

Material (material): ISO 5832-1
Plattendicke (thickness): 4,0 mm
Plattenbreite (width): 11,5 mm
Lochabstand Mittel (hole spacing middle): 35,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8020-5008	126	8
8020-5009	136	9
8020-5010	150	10
8020-5012	174	12



Veterinär Platten / veterinary plates

3,5 mm Vet -Abstützplatte Platte

3,5 mm Butress Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,0 mm
Plattenbreite (width):	9,1 mm
Kopfbreite (head width):	20,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8020-0404	60	4
8020-0406	80	6



3,5 mm Vet -Abstützplatte Platte

3,5 mm Butress Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,2 mm
Plattenbreite (width):	12,0 mm
Kopfbreite (head width):	21,7 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8020-0508	106	8
8020-0510	130	10



3,5 mm TPLO-Platte

3,5 mm TPLO Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	3,0 mm
Plattenbreite (width):	12,0 mm
Kopfbreite (head width):	24,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
8020-0603	58	3
8020-0604	79	4



2,7 mm TPLO-Platte

2,7 mm TPLO Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	2,7 mm
Plattenbreite (width):	9,0 mm
Kopfbreite (head width):	16,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes	Ausrichtung Direction
8032-0103	43	3	rechts / right
8032-0203	43	3	links / left





3,5 mm TPLO-Platte

3,5 mm TPLO Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	4,0 mm
Plattenbreite (width):	11,0 mm
Kopfbreite (head width):	18,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes	Ausrichtung Direction
8034-0103	54	3	rechts / right
8034-0203	54	3	links / left



3,5 mm TPLO-Platte

3,5 mm TPLO Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	4,0 mm
Plattenbreite (width):	11,0 mm
Kopfbreite (head width):	22,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes	Ausrichtung Direction
8035-0103	63	3	rechts / right
8035-0203	63	3	links / left




4,5 mm TPLO-Platte

4,5 mm TPLO Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	4,0 mm
Plattenbreite (width):	11,0 mm
Kopfbreite (head width):	25,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes	Ausrichtung Direction
8036-0104	78	4	rechts / right
8036-0204	78	4	links / left



1,5 / 2,0 mm Arthrodesese Platte

1,5 / 2,0 mm Arthrodesis Plate

Material (material):	ISO 5832-1
Plattendicke (thickness):	1,2 / 2,0 mm
Plattenbreite (width):	5,0 / 6,3 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401168	53	9



Veterinär Platten / veterinary plates

2,0 / 2,0 mm Arthrodesese Platte

2,0 / 2,0 mm Arthrodesis Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 1,2 / 2,0 mm
 Plattenbreite (width): 5,0 / 6,3 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401169	57	9



2,0 / 2,7 mm Arthrodesese Platte

2,0 / 2,7 mm Arthrodesis Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 1,3 / 2,0 mm
 Plattenbreite (width): 5,0 / 8,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401170	75	9



2,0 / 2,7 mm Arthrodesese Platte

2,0 / 2,7 mm Arthrodesis Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 1,3 / 2,0 mm
 Plattenbreite (width): 5,0 / 8,3 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401171	75	9



2,7 / 3,5 mm Arthrodesese Platte

2,7 / 3,5 mm Arthrodesis Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 2,0 / 3,2 mm
 Plattenbreite (width): 8,0 / 10,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401172	119	9



2,7 / 3,5 mm Arthrodesese Platte

2,7 / 3,5 mm Arthrodesis Plate

Material (material): ISO 5832-1
 Plattendicke (thickness): 2,0 / 3,4 mm
 Plattenbreite (width): 8,0 / 10,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401173	119	9





2,7 / 3,5 mm; Arthrodesese Platte

2,7 / 3,5 mm; Arthrodesis Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 2,0 / 3,4 mm
Plattenbreite (width): 8,0 / 10,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401174	102	9



2,7 / 3,5 mm; Arthrodesese Platte

2,7 / 3,5 mm; Arthrodesis Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 2,0 / 3,4 mm
Plattenbreite (width): 8,0 / 10,0 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401175	102	9



3,5 / 3,5 mm; Arthrodesese Platte

3,5 / 3,5 mm; Arthrodesis Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 2,5 / 4,0 mm
Plattenbreite (width): 8,2 / 11,5 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401176	102	11



3,5 / 3,5 mm Arthrodesese Platte

3,5 / 3,5 mm Arthrodesis Plate

Material (material): ISO 5832-1
Plattendicke (thickness): 2,4 / 4,0 mm
Plattenbreite (width): 8,5 / 11,5 mm

Art.Nr. Item no.	Länge mm length mm	Löcher Holes
9107401177	104	10



1,5 mm Kortikalis Schraube, Innensechskant

1,5 mm Cortical Screw, hexagonal recess

Material (material): ISO 5832-1 /-3
 Gewindedurchmesser (thread diameter): 1,5 mm
 Kerndurchmesser (core diameter) : 1,0 mm
 Gewindesteigung (thread pitch): 0,5 mm

Art.Nr. Item no. standard ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 neutral	Länge mm length mm
6060-0006	6060-0506	6160-0506	6
6060-0007	6060-0507		7
6060-0008	6060-0508	6160-0508	8
6060-0009	6060-0509		9
6060-0010	6060-0510	6160-0510	10
6060-0011	6060-0511		11
6060-0012	6060-0512	6160-0512	12
6060-0013	6060-0513		13
6060-0014	6060-0514	6160-0514	14
6060-0016	6060-0516	6160-0516	16
6060-0018	6060-0518	6160-0518	18
6060-0020	6060-0520	6160-0520	20
6060-0022	6060-0522	6160-0522	22
6060-0024	6060-0524	6160-0524	24



1,5 mm Kortikalis Schraube, Kreuzschlitz

1,5 mm Cortical Screw, cruciform recess

Material (material): ISO 5832-1
 Gewindedurchmesser (thread diameter): 1,5 mm
 Kerndurchmesser (core diameter) : 1,0 mm
 Gewindesteigung (thread pitch): 0,5 mm

Art.Nr. Item no. self cutting	Länge mm length mm
6060-2506	6
6060-2507	7
6060-2508	8
6060-2509	9
6060-2510	10
6060-2511	11
6060-2512	12
6060-2513	13
6060-2514	14
6060-2516	16
6060-2518	18
6060-2520	20
6060-2522	22





2,0 mm Kortikalis Schraube, Innensechskant

2,0 mm Cortical Screw, hexagonal recess

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 2,0 mm
Kerndurchmesser (core diameter) : 1,3 mm
Gewindesteigung (thread pitch): 0,6 mm

Art.Nr. Item no. standard ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 neutral	Länge mm length mm
6061-0006	6061-0506	6161-0506	6
6061-0008	6061-0508	6161-0508	8
6061-0010	6061-0510	6161-0510	10
6061-0012	6061-0512	6161-0512	12
6061-0014	6061-0514	6161-0514	14
6061-0016	6061-0516	6161-0516	16
6061-0018	6061-0518	6161-0518	18
6061-0020	6061-0520	6161-0520	20
6061-0022	6061-0522	6161-0522	22
6061-0024	6061-0524	6161-0524	24
6061-0026	6061-0526	6161-0526	26
6061-0028	6061-0528	6161-0528	28
6061-0030	6061-0530	6161-0530	30
6061-0032	6061-0532	6161-0532	32
6061-0034	6061-0534		34
6061-0036	6061-0536		36
6061-0038	6061-0538		38





Kortikalis Schrauben / cortical screws

2,0 mm Kortikalis Schraube, Kreuzschlitz

2,0 mm Cortical Screw, cruciform recess

Material (material): ISO 5832-1
Gewindedurchmesser (thread diameter): 2,0 mm
Kerndurchmesser (core diameter) : 1,3 mm
Gewindesteigung (thread pitch): 0,6 mm

Art.Nr. Item no. self cutting	Länge mm length mm
6061-2506	6
6061-2508	8
6061-2510	10
6061-2512	12
6061-2514	14
6061-2516	16
6061-2518	18
6061-2520	20
6061-2522	22
6061-2524	24
6061-2526	26
6061-2528	28
6061-2530	30
6061-2532	32
6061-2534	34
6061-2536	36





2,4 mm Kortikalis Schraube, Innensechskant

2,4 mm Cortical Screw, hexagonal recess

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 2,4 mm
Kerndurchmesser (core diameter) : 1,7 mm
Gewindesteigung (thread pitch): 1,0 mm

Art.Nr. Item no. self cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1 blue	Länge mm length mm
6061-5506	6161-1506	6
6061-5507		7
6061-5508	6161-1508	8
6061-5509		9
6061-5510	6161-1510	10
6061-5511		11
6061-5512	6161-1512	12
6061-5513		13
6061-5514	6161-1514	14
6061-5516	6161-1516	16
6061-5518	6161-1518	18
6061-5520	6161-1520	20
6061-5522	6161-1522	22
6061-5524	6161-1524	24
6061-5526	6161-1526	26
6061-5528	6161-1528	28
6061-5530	6161-1530	30
6061-5532	6161-1532	32
6061-5534	6161-1534	34
6061-5536	6161-1536	36
6061-5538	6161-1538	38
6061-5540	6161-1540	40



Kortikalis Schrauben / cortical screws

2,7 mm Kortikalis Schraube, Innensechskant

2,7 mm Cortical Screw, hexagonal recess

Material (material): ISO 5832-1 / -3
 Gewindedurchmesser (thread diameter): 2,7 mm
 Kerndurchmesser (core diameter) : 1,9 mm
 Gewindesteigung (thread pitch): 1,0 mm

Art.Nr. Item no. standard ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1 gold	Länge mm length mm
6062-0006	6062-0506		6
6062-0008	6062-0508	6162-0508	8
6062-0010	6062-0510	6162-0510	10
6062-0012	6062-0512	6162-0512	12
6062-0014	6062-0514	6162-0514	14
6062-0016	6062-0516	6162-0516	16
6062-0018	6062-0518	6162-0518	18
6062-0020	6062-0520	6162-0520	20
6062-0022	6062-0522	6162-0522	22
6062-0024	6062-0524	6162-0524	24
6062-0026	6062-0526	6162-0526	26
6062-0028	6062-0528	6162-0528	28
6062-0030	6062-0530	6162-0530	30
6062-0032	6062-0532	6162-0532	32
6062-0034	6062-0534	6162-0534	34
6062-0036	6062-0536	6162-0536	36
6062-0038	6062-0538		38
6062-0040	6062-0540		40





3,5 mm Kortikalis Schraube, Innensechskant

3,5 mm Cortical Screw, hexagonal recess

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 3,5 mm
Kerndurchmesser (core diameter) : 2,4 mm
Gewindesteigung (thread pitch): 1,25 mm

Art.Nr. Item no. standard ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1	Art.Nr. Item no. standard ISO 5832-3 gold	Art.Nr. Item no. self cutting ISO 5832-3 gold	Länge mm length mm
6063-0010	6063-0510	6163-0010	6163-0510	10
6063-0012	6063-0512	6163-0012	6163-0512	12
6063-0014	6063-0514	6163-0014	6163-0514	14
6063-0016	6063-0516	6163-0016	6163-0516	16
6063-0018	6063-0518	6163-0018	6163-0518	18
6063-0020	6063-0520	6163-0020	6163-0520	20
6063-0022	6063-0522	6163-0022	6163-0522	22
6063-0024	6063-0524	6163-0024	6163-0524	24
6063-0026	6063-0526	6163-0026	6163-0526	26
6063-0028	6063-0528	6163-0028	6163-0528	28
6063-0030	6063-0530	6163-0030	6163-0530	30
6063-0032	6063-0532	6163-0032	6163-0532	32
6063-0034	6063-0534	6163-0034	6163-0534	34
6063-0036	6063-0536	6163-0036	6163-0536	36
6063-0038	6063-0538	6163-0038	6163-0538	38
6063-0040	6063-0540	6163-0040	6163-0540	40
6063-0042	6063-0542			42
6063-0044	6063-0544			44
6063-0045	6063-0545	6163-0045	6163-0545	45
6063-0546	6063-0546			46
6063-0548	6063-0548		6163-0548	48
6063-0050	6063-0550	6163-0050	6163-0550	50
6063-0055	6063-0555		6163-0555	55
6063-0060	6063-0560		6163-0560	60
6063-0065	6063-0565		6163-0565	65
6063-0070	6063-0570		6163-0570	70
6063-0075	6063-0575		6163-0575	75
6063-0080	6063-0580		6163-0580	80
6063-0085	6063-0585		6163-0585	85
6063-0090	6063-0590		6163-0590	90
6063-0095	6063-0595		6163-0595	95
6063-0100	6063-0600		6163-0600	100
6063-0105	6063-0605		6163-0605	105
6063-0110	6063-0610		6163-0610	110
6063-0115	6063-0615		6163-0615	115
6063-0120	6063-0620		6163-0620	120



Kortikalis Schrauben / cortical screws

4,5 mm Kortikalis Schraube, Innensechskant

4,5 mm Cortical Screw, hexagonal recess

Material (material): ISO 5832-1 / -3
 Gewindedurchmesser (thread diameter): 4,5 mm
 Kerndurchmesser (core diameter) : 3,0 mm
 Gewindesteigung (thread pitch): 1,75 mm

Art.Nr. Item no. standard ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 gold	Länge mm length mm
6064-0014	6064-0514	6164-0514	14
6064-0016	6064-0516	6164-0516	16
6064-0018	6064-0518	6164-0518	18
6064-0020	6064-0520	6164-0520	20
6064-0022	6064-0522	6164-0522	22
6064-0024	6064-0524	6164-0524	24
6064-0026	6064-0526	6164-0526	26
6064-0028	6064-0528	6164-0528	28
6064-0030	6064-0530	6164-0530	30
6064-0032	6064-0532	6164-0532	32
6064-0034	6064-0534	6164-0534	34
6064-0036	6064-0536	6164-0536	36
6064-0038	6064-0538	6164-0538	38
6064-0040	6064-0540	6164-0540	40
6064-0042	6064-0542	6164-0542	42
6064-0044	6064-0544	6164-0544	44
6064-0045	6064-0545	6164-0545	45
6064-0046	6064-0546	6164-0546	46
6064-0048	6064-0548	6164-0548	48
6064-0050	6064-0550	6164-0550	50
6064-0052	6064-0552	6164-0552	52
6064-0054	6064-0554	6164-0554	54
6064-0055	6064-0555	6164-0555	55
6064-0056	6064-0556	6164-0556	56
6064-0058	6064-0558	6164-0558	58
6064-0060	6064-0560	6164-0560	60
6064-0062	6064-0562	6164-0562	62
6064-0064	6064-0564	6164-0564	64
6064-0065	6064-0565	6164-0565	65
6064-0066	6064-0566	6164-0566	66
6064-0068	6064-0568	6164-0568	68
6064-0070	6064-0570	6164-0570	70
6064-0072	6064-0572	6164-0572	72
6064-0074	6064-0574	6164-0574	74
6064-0075	6064-0575	6164-0575	75
6064-0076	6064-0576	6164-0576	76
6064-0078	6064-0578	6164-0578	78
6064-0080	6064-0580	6164-0580	80
6064-0085	6064-0585	6164-0585	85
6064-0090	6064-0590	6164-0590	90
6064-0095	6064-0595	6164-0595	95





Kortikalis Schrauben / cortical screws

6064-0100	6064-0600	6164-0600	100
6064-0105	6064-0605	6164-0605	105
6064-0110	6064-0610	6164-0610	110
6064-0115	6064-0615	6164-0562	115
6064-0120	6064-0620	6164-0564	120
6064-0125	6064-0625	6164-0565	125
6064-0130	6064-0630	6164-0566	130
6064-0135	6064-0635	6164-0568	135
6064-0140	6064-0640	6164-0570	140

2,7 mm Verriegelungsschraube, Torx

2,7 mm Locking screw, stardrive

Material (material): ISO 5832-1 / -3
 Gewindedurchmesser (thread diameter): 2,7 mm
 Kerndurchmesser (core diameter) : 2,0 mm
 Gewindesteigung (thread pitch): 0,8 mm
 Kopfgewinde (head thread): M5

Art.Nr. Item no. self-cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 purple	Länge mm length mm
4082-0510	4182-0510	10
4082-0512	4182-0512	12
4082-0514	4182-0514	14
4082-0516	4182-0516	16
4082-0518	4182-0518	18
4082-0520	4182-0520	20
4082-0522	4182-0522	22
4082-0524	4182-0524	24
4082-0526	4182-0526	26
4082-0528	4182-0528	28
4082-0530	4182-0530	30
4082-0532	4182-0532	32
4082-0535	4182-0535	35
4082-0538	4182-0538	38
4082-0540	4182-0540	40
4082-0542	4182-0542	42
4082-0545	4182-0545	45
4082-0546	4182-0546	46
4082-0548	4182-0548	48
4082-0550	4182-0550	50
4082-0552	4182-0552	52
4082-0555	4182-0555	55
4082-0560	4182-0560	60
4082-0565	4182-0565	65
4082-0570	4182-0570	70





2,7 mm Verriegelungsschraube, Torx, polyaxial

2,7 mm Locking screw, stardrive, polyaxial

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 2,7 mm
Kerndurchmesser (core diameter) : 2,1 mm
Gewindesteigung (thread pitch): 0,9 mm
Kopfgewinde (head thread): M5

Art.Nr. Item no. self-cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 grey	Länge mm length mm
4082-1510	4182-1510	10
4082-1512	4182-1512	12
4082-1514	4182-1514	14
4082-1516	4182-1516	16
4082-1518	4182-1518	18
4082-1520	4182-1520	20
4082-1522	4182-1522	22
4082-1524	4182-1524	24
4082-1526	4182-1526	26
4082-1528	4182-1528	28
4082-1530	4182-1530	30
4082-1532	4182-1532	32
4082-1534	4182-1534	34
4082-1535	4182-1535	35
4082-1536	4182-1536	36
4082-1538	4182-1538	38
4082-1540	4182-1540	40
4082-1542	4182-1542	42
4082-1544	4182-1544	44
4082-1545	4182-1545	45
4082-1546	4182-1546	46
4082-1548	4182-1548	48
4082-1550	4182-1550	50
4082-1555	4182-1555	55
4082-1560	4182-1560	60
4082-1565	4182-1565	65
4082-1570	4182-1570	70





Verriegelungsschrauben / locking screws

3,5 mm Verriegelungsschraube, Torx

3,5 mm Locking screw, stardrive

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 3,5 mm
Kerndurchmesser (core diameter) : 2,9 mm
Gewindesteigung (thread pitch): 0,8 mm
Kopfgewinde (head thread): M5

Art.Nr. Item no. self-cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 blue	Länge mm length mm
4083-0510	4183-0510	10
4083-0512	4183-0512	12
4083-0514	4183-0514	14
4083-0516	4183-0516	16
4083-0518	4183-0518	18
4083-0520	4183-0520	20
4083-0522	4183-0522	22
4083-0524	4183-0524	24
4083-0526	4183-0526	26
4083-0528	4183-0528	28
4083-0530	4183-0530	30
4083-0532	4183-0532	32
4083-0534	4183-0534	34
4083-0535	4183-0535	35
4083-0536	4183-0536	36
4083-0538	4183-0538	38
4083-0540	4183-0540	40
4083-0542	4183-0542	42
4083-0544	4183-0544	44
4083-0545	4183-0545	45
4083-0546	4183-0546	46
4083-0548	4183-0548	48
4083-0550	4183-0550	50
4083-0552	4183-0552	52
4083-0554	4183-0554	54
4083-0555	4183-0555	55
4083-0556	4183-0556	56
4083-0558	4183-0558	58
4083-0560	4183-0560	60
4083-0562	4183-0562	52
4083-0564	4183-0564	64
4083-0565	4183-0565	52
4083-0566	4183-0566	66
4083-0568	4183-0568	68
4083-0570	4183-0570	70





3,5 mm Verriegelungsschraube, Torx, polyaxial

3,5 mm Locking screw, stardrive, polyaxial

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 3,5 mm
Kerndurchmesser (core diameter) : 2,9 mm
Gewindesteigung (thread pitch): 0,9 mm
Kopfgewinde (head thread): M5

Art.Nr. Item no. self-cutting ISO 5832-1	Art.Nr. Item no. self cutting ISO 5832-3 grau	Länge mm length mm
4083-1510	4183-1510	10
4083-1512	4183-1512	12
4083-1514	4183-1514	14
4083-1516	4183-1516	16
4083-1518	4183-1518	18
4083-1520	4183-1520	20
4083-1522	4183-1522	22
4083-1524	4183-1524	24
4083-1526	4183-1526	26
4083-1528	4183-1528	28
4083-1530	4183-1530	30
4083-1532	4183-1532	32
4083-1534	4183-1534	34
4083-1535	4183-1535	35
4083-1536	4183-1536	36
4083-1538	4183-1538	38
4083-1540	4183-1540	40
4083-1542	4183-1542	42
4083-1544	4183-1544	44
4083-1545	4183-1545	45
4083-1546	4183-1546	46
4083-1548	4183-1548	48
4083-1550	4183-1550	50
4083-1552	4183-1552	52
4083-1554	4183-1554	54
4083-1555	4183-1555	55
4083-1556	4183-1556	56
4083-1558	4183-1558	58
4083-1560	4183-1560	60
4083-1562	4183-1562	52
4083-1564	4183-1564	64
4083-1565	4183-1565	52
4083-1566	4183-1566	66
4083-1568	4183-1568	68
4083-1570	4183-1570	70



3,5 mm Spongiosa Schrauben, Innensechskant

3,5 mm Cancellous screws, hexagonal recess

Material (material): ISO 5832-1 / -3
 Gewindedurchmesser (thread diameter): 3,5 mm
 Kerndurchmesser (core diameter) : 1,8 mm
 Gewindesteigung (thread pitch): 1,75 mm

Art.Nr. Item no. self-cutting ISO 5832-1 full thread	Art.Nr. Item no. self cutting ISO 5832-3 gold	Länge mm length mm
6071-0008	6071-0008	8
6071-0010	6071-0010	10
6071-0012	6071-0012	12
6071-0014	6071-0014	14
6071-0016	6071-0016	16
6071-0018	6071-0018	18
6071-0020	6071-0020	20
6071-0022	6071-0022	22
6071-0024	6071-0024	24
6071-0026	6071-0026	26
6071-0028	6071-0028	28
6071-0030	6071-0030	30
6071-0032	6071-0032	32
6071-0034	6071-0034	34
6071-0035	6071-0035	35
6071-0036	6071-0036	36
6071-0038	6071-0038	38
6071-0040	6071-0040	40
6071-0042	6071-0042	42
6071-0044	6071-0044	44
6071-0045	6071-0045	45
6071-0046	6071-0046	46
6071-0048	6071-0048	48
6071-0050	6071-0050	50
6071-0055	6071-0055	55
6071-0060	6071-0060	60





Spongiosa Schrauben / cancellous screws

4,0 mm Spongiosa Schrauben, Innensechskant

4,0 mm Cancellous screws, hexagonal recess

Material (material): ISO 5832-1 / -3
Gewindedurchmesser (thread diameter): 4,0 mm
Kerndurchmesser (core diameter) : 1,9 mm
Gewindesteigung (thread pitch): 1,75 mm

Art.Nr. Item no. self-cutting ISO 5832-1 full thread	Art.Nr. Item no. self cutting ISO 5832-1 short thread	Art.Nr. Item no. self-cutting ISO 5832-3 full thread	Art.Nr. Item no. self cutting ISO 5832-3 short thread	Länge mm length mm
6072-0010	6072-2010	6072-0010	6072-2010	10
6072-0012	6072-2012	6072-0012	6072-2012	12
6072-0014	6072-2014	6072-0014	6072-2014	14
6072-0016	6072-2016	6072-0016	6072-2016	16
6072-0018	6072-2018	6072-0018	6072-2018	18
6072-0020	6072-2020	6072-0020	6072-2020	20
6072-0022	6072-2022	6072-0022	6072-2022	22
6072-0024	6072-2024	6072-0024	6072-2024	24
6072-0026	6072-2026	6072-0026	6072-2026	26
6072-0028	6072-2028	6072-0028	6072-2028	28
6072-0030	6072-2030	6072-0030	6072-2030	30
6072-0032	6072-2032	6072-0032	6072-2032	32
6072-0034	6072-2034	6072-0034	6072-2034	34
6072-0035	6072-2035	6072-0035	6072-2035	35
6072-0036	6072-2036	6072-0036	6072-2036	36
6072-0038	6072-2038	6072-0038	6072-2038	38
6072-0040	6072-2040	6072-0040	6072-2040	40
6072-0042	6072-2042	6072-0042	6072-2042	42
6072-0044	6072-2044	6072-0044	6072-2044	44
6072-0045	6072-2045	6072-0045	6072-2045	45
6072-0046	6072-2046	6072-0046	6072-2046	46
6072-0048	6072-2048	6072-0048	6072-2048	48
6072-0050	6072-2050	6072-0050	6072-2050	50
6072-0055	6072-2055	6072-0055	6072-2055	55
6072-0060	6072-2060	6072-0060	6072-2060	60





**ORTHO
SOLUTIONS**
IMPLANTATE TECHNOLOGIE



Ortho Solutions GmbH
Eltastraße 2
D-78573 Wurmlingen

Phone +49(0)7461 - 96357 - 0
Fax +49(0)7461 - 96357 - 99
www.orthosolutions.de
info@orthosolutions.de